



**PT-660  
CENTRIFUGE MACHINE**

**PT-660,  
PT-661,  
PT-677 &  
PT-678**

Centrifuge machines are used for:

- Laboratory Research: Separating components of a liquid or solid mixture based on density differences.
- Medical Diagnostics: Analyzing blood, urine, and other bodily fluids for diagnostic purposes.
- Biochemistry: Purifying proteins, DNA, RNA, and other biomolecules.
- Industrial Processes: Separating particles or components in manufacturing and production processes.
- Environmental Testing: Analyzing soil, water, and air samples for research and regulatory purposes.
- Quality Control: Testing and analyzing materials in various industries to ensure consistency and purity.

## PT-660 CENTRIFUGE MACHINE

Centrifuge Machine Microprocessor based, digital. Body made of rust proof ABS Plastic injection moulding with metallic protection bowl. Lid Lock for safe drive. Equipped with 0-99 minute timer. 2lines of 16 character LCD display of r.p.m., R.C.F., for set and run time. Programmable speed regulator from 500 to 4400 RPM in interval of 100 RPM. Last test parameter recall and displays automatically on LCD.

(HSN Code: 84211999)

29950-00/-

### TECHNICAL PARAMETERS WITH SWING OUT HEAD 15 ML 8 TUBE

MOTOR	Universal Motor ( with carbon brushes) 1/8 H.P
DIMENSIONS OF OUTER BOX	510 x 510 x 390 mm
CURRENT	220 Volts, 50 Hz A.C

## PT-661 CENTRIFUGE MACHINE

Centrifuge Machine Square Body Acrylic Top provided with 60 minutes timer and R.P.M controller with step control.

(Capacity of Tube 15 ml x 6 Tubes)

13500-00/-

## CENTRIFUGE TUBE

CATALOGUE NO.	CAPACITY	HT. X O.D	RATE (₹)
<b>PT-680 Centrifuge Tubes : Conical bottom, plain and borosil glass (HSN: 70179090)</b>			
PT-680/05	3 ML	76 X 10 MM	25-00/-
PT-680/10	5 ML	100 X 13 MM	30-00/-
PT-680/15	15 ML	120 X 17 MM	40-00/-
PT-680/20	50 ML	125 X 28 MM	130-00/-
<b>PT-682 Graduated, Conical bottom, and borosil glass (HSN: 70179090)</b>			
PT-682/05	5 ML	100 X 13 MM	60-00/-
PT-682/10	15 ML	120 X 17 MM	75-00/-
PT-682/15	50 ML	125 X 28 MM	155-00/-
PT-685 Conical bottom with screw cap, plain, cap. 15 ml, borosil glass			62-00/-
<b>PT-689 Conical bottom, plain, borosil glass</b>			
PT-689/05	50 ML	30 X 110 MM	110-00/-
PT-689/10	100 ML	45 X 110 MM	300-00/-
<b>PT-690 Round bottom, graduated, Borosil Glass (HSN: 70179090)</b>			
PT-690/05	50ML	30 X 110 MM	140-00/-
PT-690/10	100 ML	45 X 110 MM	385-00/-