

PT-435 KJELDHAL DIGESTION UNIT

(HSN Code: 84198990)

Temperature Range	AMBIENT + 10°C to 400°C
Temperature Accuracy	± 5°C
Temperature Controller	Energy Regulator
Construction	MS Powder coated
No. of Heaters	500 watt x 6
No. of Regulators	6
Fume duct	SS 304 with 6 holes
Stand for fume duct	MS powder coated
Kjeldhal Flask	500 ml long neck
Power supply	220 Volts 50 Hz

CATALOGUE NO.	DESCRIPTION	RATE (₹)
PT-435/05	Unit of 3 Test with Glass Parts (500 ml flask)	16950-00/-
PT-435/10	Unit of 3 Test with 3 Regulator without glass parts	14250-00/-
PT-435/15	Unit of 6 Test with Glass Parts (500 ml flask)	24950-00/-
PT-435/20	Unit of 6 Test with 6 Regulator without glass parts	22950-00/-



PT-435
KJELDHAL DIGESTION UNIT

PT-436 GLASS PARTS KJELDHAL FLASK

(HSN-70179090)

CATALOGUE NO.	DESCRIPTION	RATE (₹)
PT-436	GLASS PARTS for above Kjeldhal flasks for digestion unit borosilicate glass.	
PT-436/05	300 ml	495-00/-
PT-436/15	500 ml	725-00/-
PT-436/15	800 ml	895-00/-
PT-436/30	Fume Duct made of Stainless Steel (3 test)	2750-00/-
PT-436/35	Fume Duct made of Stainless Steel (6 test)	3850-00/-

PT-437 KJELDHAL DRY BLOCK DIGESTION

KJELDHAL DRY BLOCK DIGESTION (40 TUBES) The dry block kjeldhal is an economical, compact, constant temperature heating system suitable for 40 digestion tubes of 75ml. Dry blocks are particularly suitable in the laboratory for digestion, incubation, boiling & variety of industrial applications. The system comprises of a base heated dry block of metal alloy with energy regulator to control the temperature up to 400°C with pyrometer for viewing the temperature Heaters have been fixed at the beneath of the dry block duly embossed. A aluminum rack is provided to protect the test tubes. Without glass digestion tubes Power 2K.W., Single Phase, 230 Volts 50 Hz.

(HSN-84198990)

32500-00/-

PT-437/20

Glass digestion tubes Cap: 75 ml Borosilicate Glass (25cm × 24mm).

225-00/-

PT-440 KJELDHAL DISTILLATION UNIT

(HSN-84198990)

Temperature Range	AMBIENT + 10°C to 400°C
Temperature Accuracy	± 5°C
Temperature Controller	Energy Regulator
Construction	MS Powder coated
No. of Heaters	500 watt x 6
No. of Regulators	6
Water Circulation Assembly	Made of Brass
Stand	MS pipe powder coated
Glass Assembly	Kjeldhal assembly Borosilicate make 500 ml
Power supply	220 Volts 50 Hz

CATALOGUE NO.	DESCRIPTION	RATE (₹)
PT-440/05	Unit of 3 test, 3 regulators and 3 glass assemblies	24950-00/-
PT-440/10	Unit of 3 test, 3 regulators without glass assemblies	17950-00/-
PT-440/15	Unit of 6 test, 6 regulator and 6 glass assemblies	39950-00/-
PT-440/20	Unit of 6 test, 6 regulators without glass assemblies	24950-00/-



PT-440
KJELDHAL DIGESTION UNIT

PT-442 DISTILLATION ASSEMBLY

DISTILLATION ASSEMBLY FOR PT-440 Kjeldhal Distillation Assmebly with interchangeable glass joints size B-24. Consisting of kjeldhal flask cap. 500ml., splash head, receiver adapter and receiver flask.

(HSN Code: 84198990)

3450-00/-

PT-442 SPARE HEATER PLATE WITH HEATER

(HSN Code: 84198990)

900-00/-

PT-445 KJELDHAL DIGESTION & DISTILLATION UNIT

One single mild steel tubular stand is provided with two heating shelves on above the other on lower portion. A fume duct and a heater box are supplied for digestion part. The heaters have individual energy regulator for temperature control. Another heater box with condensor rack made of brass rods and tubes are given for distillation part. The whole unit is electrically operated on 240 Volts AC mains, complete with cord and plug. (HSN Code: 84198990)

CATALOGUE NO.	DESCRIPTION	RATE (₹)
PT-445/05	Unit of 6 test, 6 regulator and 6 glass assemblies	64500-00/-
PT-445/10	Unit of 6 test, 6 regulators without glass assemblies	42950-00/-
PT-445/15	Unit of 3 test, 3 regulators and 3 glass assemblies	39950-00/-
PT-445/20	Unit of 3 test, 3 regulators without glass assemblies	29950-00/-

PT-448 COD DIGESTION APPRATUS

This unit is supplied with fifteen sets of glass reaction vessels 200ml and air condensers with digital temperature controller cum display. Temperature range up to 250°C. It is provided with timer up to 2 hrs to work on 220 volts AC 15 Hz. Phase-1

(HSN Code: 84198990)

57950-00/-

PT-448/05
Spare glass part set for COD distillationn apratus



PT-448
COD DISTILLATION APPRATUS

