

PT-1520
CLINICAL MICROSCOPE

(HSN Code: 90118000)

MODEL	PT- 1520
Eye Pieces	H 10x and H 15x
Objectives	Achromatic 10x and 45x
Magnification	100x to 675 x
Rate	7500-00/-

PT-1550
PATHOLOGICAL/ MEDICAL MICROSCOPE

PATHOLOGICAL/MEDICAL MICROSCOPE With attachable mechanical stage for X and Y movements upto 50 and 75mm. Sub-stage condenser N.A 1.3 and iris diaphragm moveable rack and pinion motion.

(HSN-90118000)

Eye Pieces	H 5x, H10x or H 15x
Objectives	Achromatic 10x, 45x and 100x (SL) Oil immersion
Magnification	50x to 1000 x
Rate	11950-00/-

PT-1580
BINOCULAR RESEARCH MICROSCOPE

BINOCULAR RESEARCH MICROSCOPE Advanced, with attachable mechanical stage, binocular head rotatable to 360° with built in light & transformer.

(HSN-90118000)

Eye Pieces	H 5x or H 15x & H10x (W/F paired)
Objectives	Achromatic 5x, 10x, 45x and 100x (SL) Oil immersion
Magnification	50x to 1500 x
Rate	19950-00/-



PT-1800
MICROTOME ROTATORY

MICROTOME ROTARY As per Erma Japan, complete with razor etc.

(HSN-90279010)

35900-00/-

PT-1810
SENIOR ROTATORY MICROTOME

SENIOR ROTARY MICROTOME As per Spencer 820, complete with razor etc.

(HSN-90279010)

42950-00/-

MODEL	DESCRIPTION	RATE
PT-1810/05	Automatic Razor Sharpner Spencer Type	42950-00/-
PT-1810/10	Microtome knife made of carbon steel with back handle, size 120mm	2950-00/-
PT-1810/15	Microtome knife made of carbon steel with back handle, size 160 mm	3950-00/-



PT-1810
SENIOR ROTATORY MICROTOME

