

## PT-210 PULVERIZER

LABORATORY PULVERIZER

TECHNICAL SPECIFICATIONS

165000/-

Motor Capacity	3 HP
Power supply	3 Phase, 50 Cycles, 440 Volts AC
Disc Dia	175 mm
Feed Size	8 mm
Finished Product	100-150 Mesh

## PT-211 JAW CRUSHER

MECHANICAL JAW CRUSHER

MODEL NO.	PT-211/05	PT-211/10	PT-211/15
Rate	255000/-	145000/-	120000/-
Motor unit	7.5 HP 3 Phase 440 VOLTS Make - Crompton	5 HP 3 Phase 440 VOLTS Make - Crompton	3 HP 3 Phase 440 VOLTS Make - Crompton
Jaw crusher size	10x12 inch	6x8 inch	4x6 inch
Input Size	100 to 150 mm	50 to 100 mm	50 mm
Output size	10 to 20 mm	5 to 20 mm	2 to 15 mm
Capacity (per hour)	1.5 Ton	500 to 600 kg	100 to 200 kg



PT-210  
PULVERIZER



PT-211  
JAW CRUSHER

## PT-215 STUFFING BOX GLAND

All teflon/P.T.F.E. without pulley.  
(HSN-39269099)

CAT NO.	ITEM DESCRIPTION	RATE (₹)	
	SIZE OF CONE	TO FIT STIRRER SHAFT	
PT-215/10	B-19	6 mm	675 -00/-
PT-215/15	B-19	8 mm	675-00/-
PT-215/20	B-24	6 mm	895-00/-
PT-215/25	B-24	8 mm	895-00/-
PT-215/30	B-29	8 mm	1100-00/-
PT-215/35	B-29	10 mm	1100-00/-
PT-215/40	B-34	8 mm	1300-00/-
PT-215/45	B-34	10 mm	1300-00/-

## PT-220 TEFLON (P.T.F.E.) BLADE STIRRERS

TEFLON (P.T.F.E.) BLADE STIRRERS: The teflon blades can be removed and are mutually interchangeable with the glass shafts, these factors simplify assembly, cleaning and also reduce breakages. These are designed to sweep close to vessel walls and will not cause scratching.  
(HSN-70179090)

CAT NO.	ITEM DESCRIPTION				RATE (₹)
	APPROX. SHAFT DIA. (IN MM)	SHAFT LENGTH (IN MM)	SMALLEST NECK	BLADE THROUGH SOCKET	
PT-220/05	6	450	B-19	52	500-00/-
PT-220/10	8	450	B-19	52	500-00/-
PT-220/15	6	450	B-24	76	650-00/-
PT-220/20	8	450	B-24	76	650-00/-
PT-220/25	8	600	B-24	76	700-00/-
PT-220/30	10	450	B-24	76	700-00/-
PT-220/35	8	450	B-29	100	795-00/-
PT-220/40	8	600	B-29	100	795-00/-
PT-220/45	10	600	B-34	100	850-00/-
PT-220/50	10	600	B-34	125	1100-00/-
PT-220/55	10	600	B-50	150	1250-00/-

### PT-221 RUBBER BUSH

For Stirring (6 or 8 mm.).  
(HSN-70179090)

195-00/-

### PT-222

TEFLON (P.T.F.E.) BLADED STIRRERS  
All Teflon Bladed Stirrer 10mm x 450mm with 75 mm Blade.  
(HSN-70179090)

1450-00/-

