



**PT-6500
LABORATORY RUBBER TUBING**

PT-6500 & PT-6510

Rubber tubes have various applications, including:

- **Laboratory Use:** Transporting liquids or gases between apparatus in scientific experiments.
- **Medical Devices:** Connecting various medical instruments and equipment in healthcare settings.
- **Automotive Industry:** Used in engines and fuel systems for fluid transfer.
- **Industrial Machinery:** Transferring fluids, chemicals, and gases in manufacturing processes.
- **Household Applications:** Used in plumbing and irrigation systems for water and fluid transport.
- **Construction:** Employed in construction projects for pneumatic systems and water transport.

PT-6500 LABORATORY RUBBER TUBING

It is extremely soft and conforms to the Govt. specifications. In addition to softness it is resistant to dilute Acid, Alkalies, Ammonia etc.

(HSN-40091100)

MODEL	APPROX. BORE (MM)	WALL THICKNESS (MM)	PRICE/10 METER
PT-6500/05	5	1.5	350-00/-
PT-6500/10	6	1.5	450-00/-
PT-6500/15	7	1.5	500-00/-
PT-6500/20	8	1.5	595-00/-
PT-6500/25	9	2	695-00/-
PT-6500/30	10	2	850-00/-
PT-6500/35	12	2.5	975-00/-

PT-6510 PRESSURE RUBBER TUBING

Superior quality (minimum packing 10 meter)


(HSN-40091100)

MODEL	APPROX. BORE (MM)	WALL THICKNESS (MM)	PRICE/10 METER
PT-6510/05	6	3	750-00/-
PT-6510/10	6	5	975-00/-
PT-6510/15	8	3	795-00/-
PT-6510/20	8	5	1100-00/-
PT-6510/25	9	5	1250-00/-
PT-6510/30	9	3	900-00/-
PT-6510/35	10	6	1995-00/-
PT-6510/40	12	6	2150-00/-

PT-6520
SILICONE TUBING
(MINIMUM PACKING 5 METER)

Silicone tubing is suitable upto -50°C to 250°C this is autoclavable and odourless. So it is suitable for transfer of pharmaceutical liquids beverages and other food preparations. It is not irritant to skin. Due to highly elasticity it is also useful with penestatic pumps.
(HSN-40091100)

MODEL	INNER DIA. (MM)	OUTER DIA. (MM)	WALL THICKNESS (MM)	PRICE/METER
PT-6520/05	2	4	1.0	65-00/-
PT-6520/10	3	5	1.0	95-00/-
PT-6520/15	4	6	1.0	120-00/-
PT-6520/20	5	7	1.0	150-00/-
PT-6520/25	4	7	1.5	150-00/-
PT-6520/30	6	8	1.0	195-00/-
PT-6520/35	6	9	1.5	220-00/-
PT-6520/40	7	10	1.5	275-00/-
PT-6520/45	6	10	2.0	275-00/-
PT-6520/50	8	11	1.5	295-00/-
PT-6520/55	7	11	2.0	295-00/-
PT-6520/60	8	12	2.0	335-00/-
PT-6520/65	10	14	2.0	500-00/-
PT-6520/70	12	16	2.0	550-00/-
PT-6520/75	14	20	3	750-00/-
PT-6520/80	16	22	3	995-00/-



PT-6520
SILICONE TUBING



PT-6530
RUBBER CORKS
(MINIMUM 100 PIECES)
(HSN-40169990)

MODEL	TOTAL DIA. (MM) ±0.2	BOTTOM DIA. (MM) ±0.2	HEIGHT (MM) ±0.2	NO.	RATE PER 100 NOS.
PT-6530/05	9	6	16	000	250-00/-
PT-6530/10	13.0	9.0	15.5	0	295-00/-
PT-6530/15	14.0	11.0	18.5	1	350-00/-
PT-6530/20	16.5	13.0	20.0	2	395-00/-
PT-6530/25	19.0	14.0	21.0	3	495-00/-
PT-6530/30	21.0	15.0	22.0	4	550-00/-
PT-6530/35	22.5	17.5	23.0	5	600-00/-
PT-6530/40	24.5	20.0	24.0	6	795-00/-
PT-6530/45	26.0	20.0	27.5	7	890-00/-
PT-6530/50	26.5	21.0	27.5	8	995-00/-
PT-6530/55	28.0	21.5	31.0	9	1275-00/-
PT-6530/60	30.0	26.0	33.0	10	1595-00/-
PT-6530/65	33.5	30.0	33.5	11	1795-00/-
PT-6530/70	34.5	30.0	34.5	12	1995-00/-
PT-6530/75	38.5	31.0	37.0	13	2800-00/-
PT-6530/80	41.5	32.0	42.5	14	3150-00/-
PT-6530/85	43.5	34.0	40.5	15	3995-00/-
PT-6530/90	45.0	35.0	36.5	16	4250-00/-
PT-6530/95	49.5	40.0	40.5	17	4950-00/-
PT-6530/100	53.5	48.0	33.5	18	5750-00/-
PT-6530/105	58.0	44.0	45.0	19	6950-00/-
PT-6530/110	60.0	51.0	56.5	20	9500-00/-



PT-6540
SILICONE RUBBER CORKS (TRANSLUCENT)
(HSN-40169980)

MODEL	NO.	BOTTOM DIA. (MM) ±0.2	TOP DIA. (MM) ±0.2	HEIGHT (MM) ±0.2	RATE PER PIECE
PT-6540/05	9	9	11.5	20	19-00/-
PT-6540/10	11	11	14	24	25-00/-
PT-6540/15	13	14	16	24	30-00/-
PT-6540/20	15	15	18	24	40-00/-
PT-6540/25	17	17	20	26	50-00/-
PT-6540/30	18	18	21	26	58-00/-
PT-6540/35	19	19	22.5	28	65-00/-
PT-6540/40	23	23	25.5	28	85-00/-
PT-6540/45	25	25	28.5	28	100-00/-
PT-6540/50	27	27	31	32	125-00/-
PT-6540/55	29	29	33	32	160-00/-
PT-6540/60	31	31	36	35	210-00/-
PT-6540/65	33	33	38.5	38.5	295-00/-

PT-6550
RUBBER CONES FOR BUCHNER FUNNEL
(HSN-40169990)

MODEL	SIZE SUITABLE FOR	RATE
PT-6550/05	250 ml flask	85-00/-
PT-6550/10	500 ml flask	85-00/-
PT-6550/15	1000 ml flask	165-00/-
PT-6550/20	2000 ml flask	165-00/-
PT-6550/25	5000 ml flask	230-00/-
PT-6550/30	10,000 ml flask	425-00/-



PT-6560
SUB SEAL STOPPERS

Rubber autoclaveble with town over flange and serrations for test tube.
(HSN-40169990)

MODEL	TEST TUBE O.D mm	RATE
PT-6560/05	12	26-00/-
PT-6560/10	15	29-00/-
PT-6560/15	18	45-00/-
PT-6560/20	25	55-00/-

PT-6570
PIPETTE BULB (BEST QUALITY RUBBER)
(HSN-40169980)

MODEL	DESCRIPTION	RATE
PT-6570/05	Pipette Bulb Small	50-00/-
PT-6570/10	Pipette Bulb Large	60-00/-
PT-6570/15	Pipette Filler	175-00/-

PT-6580
CHLOROFORM BELLOW WITH COTTON NET
Rs. 350-00/-

PT-6590
AUTOMIZER BULB (WITH LONG TUBE)
(HSN-40169980) Rs. 175-00/-

