

PT-585 SIEVE SHAKER (GYRATORY TYPE)

Driven by a 0.25 H.P. motor, 220 volt A.C. through a reduction gear and is suitable to carry upto 7 sieves of 6" or 8" dia. When clamped the sieve table does not rotate but is inclined from the vertical axis and the direction of inclination changes progressively in the clockwise direction. In addition to the gyratory motion of the shaker there is an upward and downward movement at a frequency of 270 r.p.m. The whole mechanism runs in oil bath and the table can be unclamped to give an additional rotary motion with timer switch.

(HSN Code: 84798200) 39950-00/-

PT-585/ 01

SIEVE SHAKER (GYRATORY TYPE) provided with 12" Dia test

(HSN Code: 84798200) 45950-00/-

PT-586 SIEVE SHAKER

(HSN-96040000)

CATALOGUE NO.	B.S.S MESH NO.	ASTM NO.	I.S.S NO.	WIDTH OF APERTURE	RATE (₹)
PT-586/ 05	4	5	480	4.75 MM	1100-00/-
PT-586/ 10	5	6	340	3.35 MM	1100-00/-
PT-586/ 15	6	7	280	2.80 MM	1100-00/-
PT-586/ 20	7	8	240	2.36 MM	1100-00/-
PT-586/ 25	8	10	200	2.00 MM	1100-00/-
PT-586/ 30	10	12	170	1.70 MM	1100-00/-
PT-586/ 35	12	14	140	1.40 MM	1100-00/-
PT-586/ 40	14	16	120	1.18 MM	1100-00/-
PT-586/ 45	16	18	100	1.00 MM	1100-00/-
PT-586/ 50	18	20	85	850 Microns	1100-00/-
PT-586/ 55	22	25	70	710 Microns	1100-00/-
PT-586/ 60	25	30	60	600 Microns	1100-00/-
PT-586/ 65	30	35	50	500 Microns	1100-00/-
PT-586/ 70	36	40	40	425 Microns	1100-00/-
PT-586/ 75	44	45	45	355 Microns	1100-00/-
PT-586/ 80	52	50	30	300 Microns	1100-00/-
PT-586/ 85	60	60	25	250 Microns	1100-00/-
PT-586/ 90	72	70	20	212 Microns	1100-00/-
PT-586/ 95	85	80	18	180 Microns	1100-00/-
PT-586/ 100	100	100	15	150 Microns	1100-00/-
PT-586/ 105	120	120	12	125 Microns	1200-00/-
PT-586/ 110	150	140	10	106 Microns	1300-00/-

PT-586/ 115	170	170	9	90 Microns	1400-00/-
PT-586/ 120	200	200	8	75 Microns	1500-00/-
PT-586/ 125	240	230	7	63 Microns	1600-00/-
PT-586/ 130	300	270	6	53 Microns	1700-00/-
PT-586/ 135	325	300	5	45 Microns	1800-00/-
PT-586/ 140	350	325	4	44 Microns	2200-00/-
PT-586/ 145	400	380	3	40 Microns	2650-00/-

PT-587 LID AND PAN SET

LID and PAN set made of brass for above

1100-00/-

PT-595 COARSE SIEVES

COARSE SIEVES 12' dia, (30 m), GI sheet frame, with steel sheet square punching.

(HSN Code: 96040000) 750-00 /-

125 mm	106 mm	100 mm	90 mm	80 mm
75 mm	63 mm	53 mm	50 mm	45 mm
40 mm	37.5 mm	31.5 mm	26.5 mm	25 mm
22.4 mm	20 mm	16 mm	13.2 mm	12.5 mm
11.2 mm	10 mm	8 mm	6.7 mm	6.3 mm
5.6 mm	4.75 mm			



PT-585
SIEVE SHAKER (GYRATORY TYPE)

