

## PT-193 LABORATORY STIRRER

(HSN-84790000)

### TECHNICAL SPECIFICATIONS

RPM	200-2500	
Motor	24 V 25 W T3000 RPM	
Speed regulator	Speed regulator in separate box	
Standard fittings	- Height adjustable Stand - Cast iron base	- SS 304 Shaft - SS 304 Blade (37 mm sweep)
Power supply	220 Volts 50Hz	

CATALOGUE NO.	RPM CONTROLLER	RATE (₹)
PT-193/05	ANALOG	13950/-
PT-193/10	DIGITAL	18950/-



PT-193  
LABORATORY STIRRER



PT-194  
TISSUE HOMOGENISERS

## PT-194 TISSUE HOMOGENISERS (Low Speed)

It is suitable for low speed homogenisation of tissue and also small volume low speed stirring applications. Low speed operation provides high torque for grinding and homogenising of animal or vegetable tissues. Provided with glass cup with teflon (P.T.F.E.) pestle having stainless steel shaft, capacity 30 ml and variable speed controlled from 200 to 2400 RPM.

(HSN-84798200)

13950-00/-

### PT-196 RUBBER BUSH

For stirring ( 6 or 8mm)

195-00/-

## PT-200 TEFLON PESTLE TISSUE HOMOGENISERS (GRINDER)

Consisting of piston type teflon pestle with serrated tip fitted to a stainless steel rod and a grinding vessel of borosilicate glass for the preparation of tissue homogenates.

(HSN-84798200)

CATALOGUE NO.	WORKING CAPACITY	RATE (₹)
PT-200/05	5 ml	725/-
PT-200/10	10 ml	975/-
PT-200/15	30 ml	1225/-
PT-200/20	55 ml	1525/-
PT-200/25	100 ml	1975/-

## PT-205 TEFLON MAGNETIC ROTORS

CATALOGUE NO.	LENGTH X DIA	RATE (₹)
PT-205	TEFLON/P.T.F.E COATED MAGNETIC ROTORS/NEEDLES	
PT-205/05	12 X 8 MM	140/-
PT-205/10	18 X 8 MM	175/-
PT-205/15	25 X 8 MM	200/-
PT-205/20	30 X 8 MM	250/-
PT-205/25	35 X 8 MM	300/-
PT-205/30	50 X 9 MM	375/-
PT-205/35	Teflon Policeman 300mm	1250/-