



PT-32 TO 35  
UV INSPECTION CABINET

## PT-32 TO 35

UV inspection cabinets have several applications, including:

- Quality Control in Manufacturing: Detecting surface defects, contamination, or foreign particles on products.
- Medical Diagnostics: Analyzing samples such as urine or blood for fluorescence under UV light.
- Geology and Mineralogy: Examining rocks and minerals for fluorescent properties.
- Document Verification: Checking the authenticity of currency, passports, and other important documents.
- Sanitation Verification: Ensuring cleanliness in hygiene-sensitive environments, such as hospitals and laboratories.

## PT-32 TO 35 UV INSPECTION CABINET

(HSN:90275090)

CATALOGUE NO.	VISIBLE LIGHT	365 nm	312 nm	SHORT UV	254 nm (with filter)	RATE (₹)
PT-32	✓	✓				8950/-
PT-33	✓	✓		✓		9950/-
PT-34	✓	✓			✓	16500/-
PT-35	✓	✓	✓		✓	23950/-
<b>SPARE UV TUBES</b>						
PT-36/05			UV TUBE 365 nm			1200/-
PT-36/10			UV TUBE 254 nm			1200/-
PT-36/15			UV TUBE 312 nm			4950/-

## PT-42 ULTRA VOILET TORCH

It is a portable device for emitting U.V. light provided with two long waves U.V. source for examining fluorescent stains or spots in gels or on T.L.C. plates. U.V. source are protected under U.V. transparent window.

(HSN-85394900)

3900-00/-

## PT-43 ULTRA VOILET PROTECTIVE GOGGLES

1450-00/-

