

PT-4122
WEIGHT BOXES

ANALYTICAL WEIGHT BOX WITH CALIBRATION CERTIFICATE BY (N.A.B.L. AUTHORIZED LABORATORY) reading of every Individual weight 1m. gm to 200 gm. (Accepted by ISO 9000 Centification agencies)

(HSN-84198990) 11950-00/-

PT-4125
ANALYTICAL WEIGHT BOX

ANALYTICAL WEIGHT BOX Made of stainless steel packed in polished wooden box with fractional weights and forcep. Class 'A', with works certificate having value of each weight. (HSN-84239010)

MODEL	DESCRIPTION	RATE
PT-4125/05	Cap. 1 mgm to 20 gms	2100-00/-
PT-4125/10	Cap. 1 mgm to 50 gms	2300-00/-
PT-4125/15	Cap. 1 mgm to 100 gms	2500-00/-
PT-4125/20	Cap. 1 mgm to 200 gms	3975-00/-



PT-4125
ANALYTICAL WEIGHT BOX



PT-4130
ANALYTICAL WEIGHT BOX

Made of brass, chrome plated packed in polished wooden box with fractional weights and forcep. Class 'A' as per I.S. specifications, with works certificate having value of each weight. (HSN-84239010)

MODEL	DESCRIPTION	RATE
PT-4130/10	Cap. 1 mgm to 50 gms	2300-00/-
PT-4130/15	Cap. 1 mgm to 100 gms	2500-00/-
PT-4130/20	Cap. 1 mgm to 200 gms	3750-00/-

PT-4140
PHYSICAL WEIGHT BOX

Made of brass, nickel plated with forceps and polished wooden case. (HSN-84239010)

MODEL	DESCRIPTION	RATE
PT-4140/05	Cap. 1 mgm to 100 gms	1500-00/-
PT-4140/10	Cap. 1 mgm to 200 gms	2100-00/-
PT-4140/15	Cap. 1 mgm to 500 gms	4450-00/-
PT-4140/20	Cap. 1 mgm to 1000 gms	7500-00/-
PT-4140/25	Cap. 1 mgm to 2000 gms	13950-00/-

PT-4150

ANALYTICAL FRACTIONAL WEIGHTS

Class 'A' in a round plastic box without forcep.

(HSN-84239010) 1300-00/-

PT-4155

ANALYTICAL FRACTIONAL WEIGHTS

Class 'B' in a round plastic case.

(HSN-84239010) 400-00/-

