



POPULAR TRADERS

Your Partner in Laboratory & Pharmaceutical Solutions



ANALYTICAL WEIGHT BOX

PT-4125 Stainless Steel · Class 'A' · 1 mg to 200 g

*Stainless steel analytical weights in a polished wooden box with fractional weights and forcep.
Class 'A', with works certificate giving the value of each weight.*



CERTIFIED & COMPLIANT



Stainless Steel

Weights

Class 'A'

Weight class

1 mg - 200 g

Capacity range

Works cert.

Value of each weight

Four capacities: 1 mg to 20 g, 50 g, 100 g or 200 g · Polished wooden box



Technical Specifications

Parameter	Specification
Material	Stainless steel
Class	Class 'A'
Case	Polished wooden box
Includes	Fractional weights and forcep
Certificate	Works certificate, each weight
Capacity	1 mg to 20 / 50 / 100 / 200 g
HSN Code	84239010

Model Options

Catalogue	Capacity
PT-4125/05	1 mg to 20 g
PT-4125/10	1 mg to 50 g
PT-4125/15	1 mg to 100 g
PT-4125/20	1 mg to 200 g

Applications

- Checking laboratory balances
- Analytical weighing in chemistry laboratories
- Teaching and QC laboratories

Why PT-4125

- Stainless steel weights.
- Class 'A' with works certificate.
- Fractional weights and forcep included.
- Four capacity options.

Standards & Quality

- Manufactured under an ISO 9001:2015 quality system.

Request a Quotation

+91-9896030622 · +91-8199900060
popular1942@gmail.com