



POPULAR TRADERS

Your Partner in Laboratory & Pharmaceutical Solutions



BULK DENSITY APPARATUS

PT-1167 A With Digital Revolution Counter · 1/6 HP Motor · 50 ml Cylinder

A simple and standard means of measuring the bulk density of powders and granules.

Provided with a digital revolution counter, 1/6 HP motor and 50 ml measuring cylinder.

50
ML



CERTIFIED & COMPLIANT



MSME
MICRO, SMALL & MEDIUM ENTERPRISES
सूक्ष्म, लघु एवं मध्यम उद्यम
OUR STRENGTH • हमारी शक्ति

Ministry of MSME, Govt. of India



GeM
Government
e Marketplace

Digital

Revolution counter

1/6 HP

Motor

50 ml

Measuring cylinder

230 V AC

Single phase

Bulk density of powders, granules and similar products · Digital revolution counter



Technical Specifications

Parameter	Specification
Revolution Counter	Digital
Motor	1/6 HP
Measuring Cylinder	50 ml
Power	AC / 230 Volts / Single Phase
HSN Code	90278090

Model Options

Catalogue	Description
PT-1167	Bulk density apparatus, manual RPM
PT-1167 A	With digital revolution counter

Applications

- Bulk density of powders and granules
- Pharmaceutical formulation and QC laboratories
- Testing of similar powdered products

Why PT-1167 A

- Simple, standard bulk density method.
- Digital revolution counter.
- 1/6 HP motor drive.
- Manual version available as PT-1167.

Standards & Quality

- Manufactured under an ISO 9001:2015 quality system.

Request a Quotation

+91-9896030622 · +91-8199900060
popular1942@gmail.com