



POPULAR TRADERS

Your Partner in Laboratory & Pharmaceutical Solutions



COD DIGESTION APPARATUS

PT-448 15 Reaction Vessels · Up to 250°C · Digital Controller · 2 Hr Timer

Supplied with fifteen sets of 200 ml glass reaction vessels and air condensers.

Digital temperature controller cum display, range up to 250°C, timer up to 2 hrs.



CERTIFIED & COMPLIANT



15 Sets

Reaction vessels

200 ml

Vessel capacity

Up to 250°C

Temperature range

2 Hrs

Timer

Air condensers · Digital temperature controller cum display · 220 V AC, single phase



Technical Specifications

Parameter	Specification
Reaction Vessels	15 sets, glass, 200 ml
Condensers	Air condensers
Temperature Range	Up to 250°C
Controller	Digital temperature controller cum display
Timer	Up to 2 hrs
Power Supply	220 Volts AC, single phase
HSN Code	84198990

Optional Accessories

- PT-448/05 Spare glass part set

Applications

- Chemical Oxygen Demand (COD) testing
- Water and wastewater analysis
- Environmental and effluent testing laboratories

Why PT-448

- Fifteen samples in one run.
- Digital controller with temperature display.
- Built-in timer up to 2 hrs.
- Spare glass part set available.

Standards & Quality

- Manufactured under an ISO 9001:2015 quality system.
- MSME (Udyam) registered, ZED certified manufacturer, listed on GeM.

Request a Quotation

+91-9896030622 · +91-8199900060
popular1942@gmail.com