



POPULAR TRADERS

Your Partner in Laboratory & Pharmaceutical Solutions



CENTRIFUGE MACHINE

PT-660 Microprocessor Based Digital · 500 to 4400 RPM · 8 x 15 ml

Microprocessor based digital centrifuge with 2 x 16 character LCD showing RPM and RCF.

Rust proof ABS body with metallic protection bowl and lid lock.



CERTIFIED & COMPLIANT



500-4400

Programmable RPM

0-99 min

Timer

LCD

RPM, RCF, run time

8 x 15 ml

Swing out head

100 RPM speed steps · Last test parameter recall · Lid lock for safe drive



Technical Specifications

Parameter	Specification
Control	Microprocessor based, digital
Speed	Programmable 500 to 4400 RPM, 100 RPM steps
Timer	0-99 minute
Display	2 lines of 16 character LCD: RPM, RCF, set and run time
Head	Swing out, 15 ml x 8 tube
Body	Rust proof ABS plastic, metallic protection bowl
Safety	Lid lock
Memory	Last test parameters recalled automatically
Motor	Universal motor (carbon brushes) 1/8 H.P.
Outer Box	510 x 510 x 390 mm
Current	220 Volts, 50 Hz A.C.
HSN Code	84211999

Model Options

Catalogue	Description
PT-660	Microprocessor based digital; 15 ml x 8 tubes
PT-661	Square body, acrylic top; 60 min timer, step RPM control; 15 ml x 6 tubes

Optional Accessories

- PT-680: Centrifuge tubes, conical, borosil glass (3, 5, 15, 50 ml)
- PT-682: Graduated centrifuge tubes, borosil glass

Applications

- Laboratory research: separating mixtures by density
- Medical diagnostics: blood and urine analysis
- Biochemistry: purifying proteins, DNA and RNA
- Environmental testing of soil and water samples
- Quality control in industry

Why PT-660

- Programmable speed with last-parameter recall.
- LCD shows RPM, RCF, set and run time.
- Rust proof ABS body and lid lock.
- 0-99 minute timer.

Standards & Quality

- Manufactured under an ISO 9001:2015 quality system.

Request a Quotation

+91-9896030622 · +91-8199900060
popular1942@gmail.com