



POPULAR TRADERS

Your Partner in Laboratory & Pharmaceutical Solutions



CRUDE FIBRE APPARATUS

PT-109 Estimation · Hot Plate Type · Six Places

Six individually regulated hot plates for crude fibre estimation in one compact unit.

Supplied complete with beakers, bulb type condensers and clamps.

6
PLATES



CERTIFIED & COMPLIANT



MSME
MICRO, SMALL & MEDIUM ENTERPRISES
सूक्ष्म, लघु एवं मध्यम उद्यम

OUR STRENGTH • हमारी शक्ति

Ministry of MSME, Govt. of India



GeM
Government
e Marketplace

6 Plates

4½" dia. hot plates

900 W

Max total power rating

Individual

Energy regulators

Complete Set

Beakers, condensers, clamps

Mild steel housing · 220/230 V mains · Spare condensers and Teflon assembly available



Technical Specifications

Parameter	Specification
Housing	Mild steel
Hot Plates	Six, 4½" dia., mounted on top of housing
Control	Individual energy regulator
Max. Power Rating	900 W AC
Power Supply	220/230 volts mains
Beakers	Six tall form spoutless beakers
Condensers	Six bulb type condensers
Clamps	For holding glass parts
HSN Code	84198990

Optional Accessories

- PT-109/05 Spare condenser for crude fibre with tall form spoutless beaker
- PT-110 Teflon assembly with four glass tubes (one set) for crude fibre

Applications

- Dietary analysis: fibre content in food and feed samples
- Feed analysis: fibre content in animal feed
- Research and development on dietary fibre
- Quality assurance for regulatory compliance on fibre content
- Environmental analysis of plant materials

Why PT-109

- Six samples run side by side with separate heat control.
- Ready to use: glassware and clamps included.
- Rugged mild steel housing.

Standards & Quality

- Manufactured under an ISO 9001:2015 certified quality system.
- Popular Traders: CE marked, MSME (Udyam) registered, ZED certified, listed on GeM.
- Serving laboratories and industry since 1942.

Request a Quotation

+91-9896030622 · +91-8199900060
popular1942@gmail.com