



POPULAR TRADERS

Your Partner in Laboratory & Pharmaceutical Solutions



DIGITAL COLONY COUNTER

PT-859 Pharmaceutical · Illuminated Grid · 4 Digit Counter

Rapid, accurate counting and recording of bacterial and mould colonies.

Pen-type probe marks the petri dish and counts at the same time, with an audio beep.

9999
COUNTS



CERTIFIED & COMPLIANT



4 Digit

Resettable counter

Audio Beep

At every count

100 mm

Magnifying lens

Cool Light

Fluorescent illumination

Pen-type probe · Wolffhügel grid glass plate · Uniform glare-free illumination



Technical Specifications

Parameter	Specification
Counter	4 digit, 9999 counts, resettable
Audio	Beep at every count
Probe	Auto marker pen probe
	Marks and counts simultaneously
Grid	Wolffhügel grid glass plate
Illumination	Cool fluorescent tube, glare-free
Magnification	100 mm dia magnifying lens
Weight	4.5 kg approx.

Optional Accessories

- PT-860 Petri platter (manual)

Applications

- Counting of bacterial and mould colonies
- Biology and zoology laboratories
- Bacteriology and pathology laboratories
- Laboratories where colony counting is a regular routine

Why PT-859

- Marking and counting in one action with the pen probe.
- Audio beep confirms every count.
- Uniform, glare-free cool-light grid.
- Magnifying lens for easy viewing.

Standards & Quality

- Manufactured under an ISO 9001:2015 quality system.
- Built for routine colony counting.

Request a Quotation

+91-9896030622 · +91-8199900060
popular1942@gmail.com