



POPULAR TRADERS

Your Partner in Laboratory & Pharmaceutical Solutions



DIGITAL VARIABLE VOLUME MICRO PIPETTES

PT-5500 0.5 to 1000 μ l · Four Volume Ranges · Stream Lined Tip Ejector

Digital variable volume micro pipettes with a stream lined tip ejector.

Volume is set by turning the plunger; two step operation.



CERTIFIED & COMPLIANT



0.5-1000 μ l

Volume span

4 Models

Variable volume

Two step

Plunger operation

Stream lined

Tip ejector

Ranges 0.5-10, 5-50, 20-200 and 200-1000 μ l · Steps from 0.1 μ l · Tips and stand available



Technical Specifications

Parameter	Specification
Type	Digital, variable volume
Volume Ranges	0.5-10 / 5-50 / 20-200 / 200-1000 µl
Volume Steps	0.1 / 1.0 / 1.0 / 10.0 µl
Volume Adjustment	By turning the plunger
Operation	Two step
Tip Ejector	Stream lined
HSN Code	90318000

Model Options

Catalogue	Capacity	Steps
PT-5500/05	0.5-10 µl	0.1 µl
PT-5500/10	5-50 µl	1.0 µl
PT-5500/15	20-200 µl	1.0 µl
PT-5500/20	200-1000 µl	10.0 µl

Optional Accessories

- Yellow tips (0-200 µl), pack of 1000 - PT-5530/05
- Blue tips (200-1000 µl), pack of 1000 - PT-5530/10
- Stand for micro pipette - PT-5530/15

Applications

- Transfer of small liquid volumes
- Sample and reagent dispensing
- Clinical, research and teaching laboratories

Why PT-5500

- Digital variable volume setting.
- Four ranges from 0.5 µl to 1000 µl.
- Fine 0.1 µl steps on the 0.5-10 µl model.
- Matching tips and pipette stand available.

Standards & Quality

- Manufactured under an ISO 9001:2015 quality system.
- CE marked, MSME registered, ZED certified and listed on the GeM portal.

Request a Quotation

+91-9896030622 · +91-8199900060
popular1942@gmail.com