



POPULAR TRADERS

Your Partner in Laboratory & Pharmaceutical Solutions



FLASK SHAKER (WRIST ACTION)

PT-169 Wrist Action Type · 4 or 8 Flasks up to 1000 ml · Stepless speed

Wrist action shaker built from heavy steel sheet with epoxy based polyester paint.

Electronic stepless speed control on rubber feet that absorb vibration.



CERTIFIED & COMPLIANT



Wrist action

Shaking type

4 or 8

Flasks

1000 ml

Max. flask capacity

Stepless

Electronic speed control

Heavy steel body · Rubber feet · Supplied with side arms, clamps and key



Technical Specifications

Parameter	Specification
Type	Flask shaker, wrist action
Flask Capacity	Up to 1000 ml
No. of Flasks	4 or 8
Construction	Heavy steel sheet, epoxy based polyester paint
Mounting	Four rubber feet to absorb vibration
Speed Control	Electronic, stepless
Supplied With	Side arms, clamps and key
Motor	F.H.P. motor
Power	220 Volts AC, 50 Hz
HSN Code	84798200

Model Options

Catalogue	Description
PT-169/05	To accept 4 flask
PT-169/10	To accept 8 flask

Applications

- Shaking and mixing of samples in flasks
- Laboratory mixing and agitation work
- General research and quality control labs

Why PT-169

- Accepts 4 or 8 flasks up to 1000 ml.
- Stepless electronic speed control.
- Rubber feet absorb vibration.
- Complete with side arms, clamps and key.

Standards & Quality

- Manufactured under an ISO 9001:2015 quality system.

Request a Quotation

+91-9896030622 · +91-8199900060
popular1942@gmail.com