



POPULAR TRADERS

Your Partner in Laboratory & Pharmaceutical Solutions



FLOCCULATOR ALUM JAR TEST APPARATUS

PT-1280 6 Stirrers · 25 to 100 RPM · FHP Geared Motor

Designed for flocculation tests in water treatment plants.

Six stirrers for 1 litre vessels, variable 25 to 100 RPM, with bottom light.



CERTIFIED & COMPLIANT



6

Stirrers

25 - 100 RPM

Variable speed

1 Litre

Vessel capacity

FHP

Geared motor

Special light arrangement at bottom · Digital RPM and timer version: PT-1285



Technical Specifications

Parameter	Specification
Stirrers	6
Vessel Capacity	1 litre each
Speed	Variable, 25 to 100 RPM
Motor	FHP geared motor
Lighting	Special light arrangement at bottom
HSN Code	90318000

Model Options

Catalogue	Description
PT-1280	Variable speed 25 - 100 RPM
PT-1285	With digital RPM and digital timer

Applications

- Water treatment: best coagulant dose
- Wastewater treatment: contaminant removal efficiency
- Industry: optimising chemical use in process water
- Environmental monitoring: water quality and treatment
- Research: new treatment methods

Why PT-1280

- Six stirrers run together.
- Variable speed, 25 to 100 RPM.
- Bottom light for viewing floc formation.
- PT-1285 adds digital RPM and digital timer.

Standards & Quality

- Manufactured under an ISO 9001:2015 quality system.

Request a Quotation

+91-9896030622 · +91-8199900060
popular1942@gmail.com