



POPULAR TRADERS

Your Partner in Laboratory & Pharmaceutical Solutions

HAND REFRACTOMETER

PT-465 Four Scale Ranges · 0 - 32 % to 58 - 92 % · ERMA, Made in Japan

A hand-held refractometer for quick percentage readings of liquid samples.

Four scale ranges from 0 - 32 % up to 58 - 92 %, ERMA make, Japan.

4
RANGES



CERTIFIED & COMPLIANT



MSME

MICRO, SMALL & MEDIUM ENTERPRISES

सूक्ष्म, लघु एवं मध्यम उद्यम

OUR STRENGTH • हमारी शक्ति

Ministry of MSME, Govt. of India



GeM

Government
e Marketplace

4

Range models

0 - 32 %

Lowest range

58 - 92 %

Highest range

ERMA

Made in Japan

Ranges: 0 - 32 % · 0 - 50 % · 28 - 62 % · 58 - 92 %



Technical Specifications

Parameter	Specification
Instrument	Hand Refractometer
Make	ERMA (Made in Japan)
Range Options	0-32 %, 0-50 %, 28-62 %, 58-92 %
Catalogue Nos.	PT-465/05, /10, /15, /20
HSN Code	90275090

Model Options

Catalogue	Range
PT-465/05	0 - 32 %
PT-465/10	0 - 50 %
PT-465/15	28 - 62 %
PT-465/20	58 - 92 %

Applications

- Quick concentration checks of liquid samples
- Sugar content checks in food and beverage work
- Quality control and teaching laboratories

Why PT-465

- Compact hand-held instrument for bench or field use.
- Four ranges to match the sample.
- ERMA make, Japan.

Standards & Quality

- Supplied by an ISO 9001:2015 certified company.

Request a Quotation

+91-9896030622 · +91-8199900060
popular1942@gmail.com