



# POPULAR TRADERS

Your Partner in Laboratory & Pharmaceutical Solutions




## IMMUNO ELECTROPHORESIS APPARATUS

**PT-63** White Perspex Migration Cell · Cold Water Circulation · Up to 8 Slides

*Walled migration cell of white perspex with cold water circulation in the outer chamber.*

*Perspex carrier studies up to eight 3" x 1" agar coated micro slides at once.*

**8**  
SLIDES



### CERTIFIED & COMPLIANT



**8 Slides**

3" x 1" micro slides

**Perspex**

White, migration cell

**Cooling**

Cold water inlets

**30x30x8**

Size (cm)

**Buffer chamber with electrodes · Clear top lid with handle · Plug, terminals and cables**



## Technical Specifications

| Parameter      | Specification                          |
|----------------|----------------------------------------|
| Migration Cell | Walled, white perspex                  |
| Outer Chamber  | Inlets & outlets for cold water        |
| Inner Chamber  | Buffer chamber fitted with electrodes  |
| Lid            | Clear transparent top lid with handle  |
| Carrier        | Perspex, micro immunological accessory |
| Slide Capacity | Up to eight 3" x 1" agar coated slides |
| Supplied With  | Plug, terminals and cables             |
| Power Pack     | Not included                           |
| Size           | 30 x 30 x 8 cm                         |
| HSN Code       | 90272000                               |

## Optional Accessories

- PT-62 Stabilized D.C. power supply (power pack)

## Applications

- Immunoelectrophoresis of serum proteins
- Antigen-antibody studies on agar slides
- Clinical and research laboratories

## Why PT-63

- Studies up to eight slides at once.
- Cold water circulation in outer chamber.
- Clear lid for viewing.
- Pairs with PT-62 power supply.

## Standards & Quality

- Manufactured under an ISO 9001:2015 quality system.
- CE marked; MSME and ZED certified; listed on GeM.

## Request a Quotation

+91-9896030622 · +91-8199900060  
popular1942@gmail.com