



POPULAR TRADERS

Your Partner in Laboratory & Pharmaceutical Solutions



LABORATORY STIRRER

PT-193 200 to 2500 RPM · SS 304 Shaft & Blade · Analog / Digital

Overhead laboratory stirrer with height adjustable stand and cast iron base.

Speed regulator in a separate box; SS 304 shaft and blade.



CERTIFIED & COMPLIANT



200-2500
RPM range

SS 304
Shaft & blade

37 mm
Blade sweep

Analog / Digital
RPM controller

Height adjustable stand · Cast iron base · Separate speed regulator box



Technical Specifications

Parameter	Specification
RPM	200-2500
Motor	24 V, 25 W
Speed Regulator	In separate box
Stand	Height adjustable
Base	Cast iron
Shaft	SS 304
Blade	SS 304 (37 mm sweep)
Power Supply	220 Volts 50 Hz
HSN Code	84790000

Model Options

Catalogue	RPM Controller
PT-193/05	Analog
PT-193/10	Digital

Applications

- Stirring and mixing of liquids and solutions
- Sample preparation
- Research and teaching laboratories

Why PT-193

- Wide 200-2500 RPM speed range.
- SS 304 shaft and blade.
- Stable cast iron base with adjustable stand.
- Analog or digital RPM controller.

Standards & Quality

- Manufactured under an ISO 9001:2015 quality system.

Request a Quotation

+91-9896030622 · +91-8199900060
popular1942@gmail.com