



POPULAR TRADERS

Your Partner in Laboratory & Pharmaceutical Solutions



MAGNETIC STIRRER WITH HOT PLATE

PT-185 SS 304 Plate · 300 W · 1, 2 and 5 Litre Stirring Capacity

Stirring and heating in one compact unit with an SS 304 hot plate.

Analog or digital temperature controller, MS or SS 304 cabinet.

**300
WATT**



CERTIFIED & COMPLIANT



300 W

Heating power

SS 304

Hot plate

Up to 5 L

Stirring capacity

Analog / Digital

Temp. controller

Plate sizes 150 x 150 / 190 x 190 mm · Powder coated MS or SS 304 cabinet · NABL calibration optional



Technical Specifications

Parameter	Specification
Stirring Capacity	1, 2, 5 Liters
Plate MOC	SS 304
Cabinet MOC	Powder coated MS / SS 304
Temperature Controller	Analog / Digital
Dimension (W x D)	150 x 150 / 190 x 190 mm
Wattage	300 W
Power Supply	220 Volts 50 Hz
HSN Code	84198990

Model Options

Catalogue	Plate (mm)	Power	Cap.	Controller
PT-185	150 x 150	300 W	2 L	Analog
PT-186	190 x 190	300 W	5 L	Analog
PT-187	150 x 150	300 W	2 L	Digital
PT-188	190 x 190	300 W	5 L	Digital

Optional Accessories

- Digital RPM display
- NABL calibrated controller
- Any specific dimensions
- IQ, OQ & PQ documentation

Applications

- Stirring and heating of solutions
- Sample and reagent preparation
- Research and teaching laboratories

Why PT-185

- Stirs and heats in one unit.
- SS 304 hot plate surface.
- Four models: analog or digital, two plate sizes.
- PT-187 and PT-188 have a 304 SS outer body with digital RPM display.

Standards & Quality

- Manufactured under an ISO 9001:2015 quality system.
- IQ/OQ/PQ documentation and NABL calibration available as options.

Request a Quotation

+91-9896030622 · +91-8199900060
popular1942@gmail.com