



POPULAR TRADERS

Your Partner in Laboratory & Pharmaceutical Solutions



PAPER STRIP ELECTROPHORESIS APPARATUS, HORIZONTAL

PT-60 Horizontal · 6 Paper Strips · Anode & Cathode Reservoirs

Horizontal paper strip electrophoresis apparatus with anode and cathode reservoirs.

Takes 6 paper strips, with safety cutoff switch and glass lid.

6
STRIPS



CERTIFIED & COMPLIANT



MSME
MICRO, SMALL & MEDIUM ENTERPRISES
सूक्ष्म, लघु एवं मध्यम उद्यम
OUR STRENGTH • हमारी शक्ति

Ministry of MSME, Govt. of India



GeM
Government
e Marketplace

2 x 550 ml

Buffer chamber

6 Strips

Paper strips

Cutoff

Safety switch

Glass

Lid material

Two platinum electrodes · Anode and cathode reservoirs · Glass lid



Technical Specifications

Parameter	Specification
Type	Horizontal
Buffer Chamber	2 x 550 ml
Paper Strips	6
Safety Cutoff Switch	Yes
Platinum Electrodes	Two
Lid Material	Glass
HSN Code	90272000

Model Options

Model	Type	Buffer (ml)	Strips	Cutoff	Lid
PT-60	Horizontal	2 x 550	6	Yes	Glass
PT-61	Vertical	2 x 200	3	No	Acrylic
PT-61 A	Vertical	2 x 200	6	No	Acrylic

Applications

- Protein separation by size or charge
- Clinical diagnostics
- Biochemical research
- Teaching laboratories

Why PT-60

- Horizontal format for 6 paper strips.
- Safety cutoff switch.
- Two platinum electrodes.
- Vertical PT-61 / PT-61 A also available.

Standards & Quality

- Manufactured under an ISO 9001:2015 quality system.
- CE marked; MSME and ZED certified; listed on GeM.

Request a Quotation

+91-9896030622 · +91-8199900060
popular1942@gmail.com