



POPULAR TRADERS

Your Partner in Laboratory & Pharmaceutical Solutions



ROTARY EVAPORATOR

PT-65 Quickvap Series · Diagonal Condenser · Motor Lift · up to 180 °C

Quickvap series rotary evaporator with diagonal condenser, motorised lift and LED rpm display.

Digital PID water bath with SS 304 inner body for solvent removal under reduced pressure.



180
°C BATH

CERTIFIED & COMPLIANT



20-200 rpm

DC motor, LED display

180 °C

Digital PID bath, 1250 W

Motor Lift

150 mm stroke, 6 mm/sec

Diagonal

1500 cm² condenser surface

Motorised lift · Digital PID bath up to 180 °C · SS 304 inner bath body



Technical Specifications

Parameter	Specification
Condenser Type	Diagonal
Cooling Surface	1500 cm ²
Evaporating Flask	100 mL to 2000 mL (800 mL supplied)
Receiving Flask	Max 1000 mL
Motor Principle	DC
Speed Range	20-200 rpm
RPM Display	LED
Lift	Motor, stroke 150 mm, 6 mm/sec
Voltage / Frequency	220-240 V / 50 Hz
Bath Controller	Digital (PID)
Bath Temp. Range	Ambient + 5 °C to 180 °C
Heat Output	1250 W
Temp. Uniformity	±1 °C
Bath Filling Volume (Max.)	3.5 L
Bath Body	SS 304 inner / Al. powder coated outer
Dimensions (L x W x H)	940 x 510 x 560 mm
Weight	31 kg
HSN Code	84198990

Model Options

Model	Condenser	Bath Ctrl.	Dimensions (mm)	Wt.
PT-65	Diagonal	Digital (PID)	940 x 510 x 560	31 kg
PT-66	Vertical	Digital (PID)	710 x 510 x 860	32 kg

Optional Accessories

- Evaporating flasks 100 mL to 2000 mL (800 mL flask supplied with unit, others optional)
- Receiving flasks up to 1000 mL

Applications

- Solvent removal and recovery in chemical and pharmaceutical labs
- Purification and concentration of solutions
- Sample preparation under reduced pressure
- Recrystallization of solid compounds
- Extraction of natural products from plant materials
- Concentrating flavours and fragrances (food and beverage)

Why PT-65

- Motorised lift with 150 mm stroke at 6 mm/sec.
- Digital PID bath heats from ambient + 5 °C to 180 °C, uniform to ±1 °C.
- Variable 20-200 rpm DC drive with LED rpm display.
- SS 304 inner bath body, aluminium powder-coated outer body.

Standards & Quality

- Manufactured under ISO 9001:2015 quality system.
- Rated for 5-40 °C ambient, 80 % relative humidity, IP 20.
- Supplied with an 800 mL evaporating flask.

Request a Quotation

+91-9896030622 · +91-8199900060
popular1942@gmail.com