



# POPULAR TRADERS

Your Partner in Laboratory & Pharmaceutical Solutions



## SEROLOGICAL WATER BATH

**PT-252** Rack Type · 8 to 28 L · Thermostatic or Digital

Deep-chamber water bath with racks for serological and tube-based work.  
Available with thermostatic or digital control, in MS or SS 304 outer body.



**CERTIFIED & COMPLIANT**



**8-28 L**

Chamber capacity

**2-6 Racks**

By model size

**Digital**

Or thermostatic control

**SS 304**

GMP outer body option

**Three chamber sizes · Thermostatic or digital controller · MS powder coated or SS 304 outer body**



## Technical Specifications

| Parameter              | Specification                         |
|------------------------|---------------------------------------|
| Product                | Serological Water Bath                |
| Chamber Capacity       | 8 L / 12 L / 28 L                     |
| Chamber Depth          | 200 mm                                |
| No. of Racks           | 2 / 4 / 6                             |
| Temperature Controller | Thermostatic or Digital               |
| Outer Body             | MS powder coated / 304 grade SS (GMP) |
| HSN Code               | 84198990                              |

## Model Options

| Catalogue | Chamber (mm) | Cap. | Racks |
|-----------|--------------|------|-------|
| PT-252/05 | 300×150×200  | 8 L  | 2     |
| PT-252/10 | 300×235×200  | 12 L | 4     |
| PT-252/15 | 500×300×200  | 28 L | 6     |

## Applications

- Serological and tube-based incubation
- Sample warming at controlled temperature
- General laboratory water bath work

## Why PT-252

- Deep chamber with racks for tube work.
- Choice of three chamber sizes.
- Digital controller option for better control.

## Standards & Quality

- Variants: MS body thermostatic; A model MS body digital; B model SS 304 body digital.
- Manufactured under ISO 9001:2015 quality system.

### Request a Quotation

+91-9896030622 · +91-8199900060  
popular1942@gmail.com