



POPULAR TRADERS

Your Partner in Laboratory & Pharmaceutical Solutions



STABILIZED D.C. POWER SUPPLY

PT-62 0 to 300 V · Constant Current up to 100 mA · Digital Control & Display

Coarse and fine control knobs to set any voltage from 0 to 300 V.

Digital control and display with short circuit and overload protection.

**300
VOLT**



CERTIFIED & COMPLIANT



0-300 V

Voltage range

100 mA

Constant current

Digital

Control & display

Protected

Short circuit / overload

Coarse & fine control knobs · Polarity reverse switch · Short circuit & overload protection



Technical Specifications

Parameter	Specification
Voltage	0 to 300 V, any setting
Current	Constant, up to 100 mA
Controls	Coarse & fine control knobs
Display	Digital control & display
Protection	Short circuit & overload
Polarity	Reverse switch
HSN Code	85362090

Applications

- Power source for gel electrophoresis
- Routine laboratory work
- Teaching and research laboratories

Why PT-62

- Set any voltage from 0 to 300 V.
- Constant current up to 100 mA.
- Digital control and display.
- Short circuit and overload protection.

Standards & Quality

- Manufactured under an ISO 9001:2015 quality system.
- CE marked; MSME and ZED certified; listed on GeM.

Request a Quotation

+91-9896030622 · +91-8199900060
popular1942@gmail.com