



POPULAR TRADERS

Your Partner in Laboratory & Pharmaceutical Solutions



UV INSPECTION CABINET

PT-32 to 35 Chromatography · Visible, 365 nm, 312 nm & 254 nm

Viewing cabinet for fluorescence inspection under visible and ultraviolet light.

Four models cover 365 nm, 312 nm and short UV 254 nm (with filter).

4
MODELS



CERTIFIED & COMPLIANT



Visible

Light on all models

365 nm

Long-wave UV, all models

312 nm

UV source, PT-35

254 nm

Short UV with filter

Models PT-32, PT-33, PT-34 and PT-35 · Spare UV tubes available



Technical Specifications

Parameter	Specification
Models	PT-32, PT-33, PT-34, PT-35
Visible Light	All models
365 nm UV	All models
312 nm UV	PT-35
Short UV	PT-33
254 nm UV (with filter)	PT-34, PT-35
Spare UV Tubes	365 nm, 254 nm, 312 nm
HSN Code	90275090

Model Options

Model	Visible	365 nm	312 nm	Short UV	254 nm
PT-32	Yes	Yes	-	-	-
PT-33	Yes	Yes	-	Yes	-
PT-34	Yes	Yes	-	-	Yes
PT-35	Yes	Yes	Yes	-	Yes

Optional Accessories

- PT-36/05 Spare UV tube 365 nm
- PT-36/10 Spare UV tube 254 nm
- PT-36/15 Spare UV tube 312 nm
- PT-42 Ultra Violet torch · PT-43 UV protective goggles

Applications

- Quality control in manufacturing: surface defects, contamination or foreign particles
- Medical diagnostics: fluorescence of samples such as urine or blood
- Geology and mineralogy: fluorescent properties of rocks and minerals
- Document verification: currency, passports and other documents
- Sanitation verification in hygiene-sensitive environments such as hospitals and laboratories

Why PT-32 to 35

- Visible and 365 nm UV together on every model.
- Step up to 312 nm or filtered 254 nm as your work needs.
- Replacement UV tubes stocked separately.

Standards & Quality

- Manufactured under an ISO 9001:2015 certified quality system.
- Popular Traders: CE marked, MSME (Udyam) registered, ZED certified, listed on GeM.
- Serving laboratories and industry since 1942.

Request a Quotation

+91-9896030622 · +91-8199900060
popular1942@gmail.com