



POPULAR TRADERS

Your Partner in Laboratory & Pharmaceutical Solutions



ULTRA VIOLET TORCH

PT-42 Portable UV Source · Two Long Wave UV Sources · UV Transparent Window

*A portable device emitting U.V. light from two long wave U.V. sources.
For examining fluorescent stains or spots in gels or on T.L.C. plates.*



CERTIFIED & COMPLIANT



Portable

Hand-held use

2 Sources

Long wave UV

UV Window

Protected sources

Gels / TLC

Spot examination

Two long wave U.V. sources · Protected under U.V. transparent window



Technical Specifications

Parameter	Specification
Type	Portable U.V. torch
UV Source	Two long wave U.V. sources
Protection	U.V. transparent window
Use	Fluorescent stains / spots in gels, T.L.C. plates
HSN Code	85394900

Optional Accessories

- PT-43 Ultra violet protective goggles

Applications

- Examining fluorescent spots on TLC plates
- Viewing fluorescent stains in gels
- Teaching and research laboratories

Why PT-42

- Portable for bench-side checks.
- Two long wave U.V. sources.
- Sources protected under a U.V. transparent window.
- Pair with PT-43 protective goggles.

Standards & Quality

- Manufactured under an ISO 9001:2015 quality system.
- CE marked; MSME and ZED certified; listed on GeM.

Request a Quotation

+91-9896030622 · +91-8199900060
popular1942@gmail.com