



POPULAR TRADERS

Your Partner in Laboratory & Pharmaceutical Solutions



UV TRANSILLUMINATOR

PT-64 Ultra Violet Transilluminator · Six 360 nm UV Sources

UV illumination from beneath the sample for gel electrophoresis, TLC and DNA work.

Fixed UV-transparent window with exhaust fan cooling and hinged eye-protective cover.



CERTIFIED & COMPLIANT



MSME
MICRO, SMALL & MEDIUM ENTERPRISES
सूक्ष्म, लघु एवं मध्यम उद्यम
OUR STRENGTH • हमारी शक्ति

Ministry of MSME, Govt. of India



GeM
Government
e Marketplace

6 x 360 nm

UV sources

25 x 13 cm

Window size

Exhaust Fan

Cooling

UV Cover

Hinged eye protection

Slab gel electrophoresis · TLC · Forensic examination · DNA analysis



Technical Specifications

Parameter	Specification
UV Source	Six 360 nm U.V. sources beneath window
Illumination	Ultraviolet from beneath the sample
Window	Fixed U.V. transparent sheet
Window Size	25 x 13 cm
Cooling	Exhaust fan
Eye Protection	Hinged U.V. obstruction cover over entire window
HSN Code	90182000

Applications

- Slab gel electrophoresis
- Thin layer chromatography
- Forensic examination
- D.N.A. analysis and experimentation

Why PT-64

- Six 360 nm sources give illumination across the full window.
- Exhaust fan keeps the unit cool in use.
- Hinged cover shields the eyes from UV and covers the whole window.

Standards & Quality

- Manufactured under an ISO 9001:2015 certified quality system.
- Popular Traders: CE marked, MSME (Udyam) registered, ZED certified, listed on GeM.
- Serving laboratories and industry since 1942.

Request a Quotation

+91-9896030622 · +91-8199900060
popular1942@gmail.com