



POPULAR TRADERS

Your Partner in Laboratory & Pharmaceutical Solutions



VERTICAL SLAB GEL ELECTROPHORESIS - UNIVERSAL

PT-58 Vertical Slab Gel Mechanism · 16 x 20 cm Gel · 13 Samples

Laboratory scale vertical slab mechanism for 16 x 20 cm gels.

Runs 13 samples simultaneously with a fixed platinum electrode assembly.



CERTIFIED & COMPLIANT



16 x 20 cm

Gel size

13 Samples

Simultaneous

Platinum

Fixed electrodes

1 mm

Comb & spacers

Notched & rectangular glass plates · Teflon comb 13 wells · Fixed gasket · Leveling screws



Technical Specifications

Parameter	Specification
Gel Size	16 x 20 cm
Samples	13 simultaneously
Electrode	Fixed platinum assembly (2 nos.)
Glass Plates	Notched & rectangular, 1 set
Comb	Teflon, 13 wells, 1 mm
Spacers	1 mm (3 nos.)
Clamps	With screws, 1 set
Gasket	Fixed
Also Included	Leveling screws, connecting cord
HSN Code	90272000

Optional Accessories

- PT-58/05 Glass plates set, 16 x 20 cm gel (spare)

Applications

- Protein separation by size or charge
- DNA and RNA analysis
- Biochemical research and clinical diagnostics
- Quality control in biotechnology

Why PT-58

- 13 samples in one run.
- Large 16 x 20 cm gel format.
- Fixed platinum electrodes.
- Spare glass plate set available.

Standards & Quality

- Manufactured under an ISO 9001:2015 quality system.
- CE marked; MSME and ZED certified; listed on GeM.

Request a Quotation

+91-9896030622 · +91-8199900060
popular1942@gmail.com