



POPULAR TRADERS

Your Partner in Laboratory & Pharmaceutical Solutions



VORTEX SHAKER

PT-158 Test Tube Shaker · Vibration Type · Variable Speed

Short, rapid eccentric stroke for vigorous agitation of test tubes and small flasks.

Shaking starts when the vessel is pressed against the cup. Variable speed, motor driven.



CERTIFIED & COMPLIANT



Variable

Speed

Eccentric

Short rapid stroke

Neoprene

Cup

Press Start

Vessel on cup

Vortex genie type · ON/OFF switch with indicator · Vacuum rubber feet prevent movement



Technical Specifications

Parameter	Specification
Type	Vibration type, vortex genie
Stroke	Short rapid eccentric stroke
Drive	Motor driven
Speed	Variable
Cup	Neoprene
Operation	Starts when vessel pressed on cup
Controls	ON/OFF switch with indicator
Base	Vacuum rubber fitted at bottom
Suitable For	Test tubes and small flasks
HSN Code	84798200

Applications

- Mixing samples in test tubes
- Mixing reagents and samples in diagnostic assays
- Clinical, microbiology and chemistry labs

Why PT-158

- Vigorous agitation from a short eccentric stroke.
- Press the tube on the cup to start shaking.
- Variable speed control.
- Rubber feet keep the unit steady on the bench.

Standards & Quality

- Manufactured under an ISO 9001:2015 quality system.
- MSME (Udyam) registered, ZED certified manufacturer, listed on GeM.

Request a Quotation

+91-9896030622 · +91-8199900060
popular1942@gmail.com